

Recommended MRI labeling based on the document, **Testing and Labeling Medical Devices for Safety in the Magnetic Resonance (MR) Environment, Guidance for Industry and Food and Drug Administration Staff**, Document issued on May 20, 2021.

**MRI Safety Information**



**MR Conditional**

The **VivoKey Spark 2** is MR Conditional. A patient with the **VivoKey Spark 2** may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.

<b>Name of the Device</b>	VivoKey Spark 2
<b>Nominal Values of Static Magnetic Field (T)</b>	1.5-T and 3-T
<b>Maximum Spatial Field Gradient (T/m and gauss/cm)</b>	40-T/m (4,000-gauss/cm)
<b>Type of RF Excitation</b>	Circularly Polarized (CP) (i.e., quadrature-driven)
<b>Transmit RF Coil Information</b>	There are no transmit RF coil restrictions. Accordingly, the following may be used: body transmit RF coil and all other RF coil combinations (i.e., body RF coil combined with any receive-only RF coil, transmit/receive head RF coil, transmit/receive knee RF coil, etc.)
<b>Operating Mode of MR System</b>	Normal Operating Mode
<b>Maximum Whole Body Averaged SAR</b>	2-W/kg (Normal Operating Mode)
<b>Limits on Scan Duration</b>	Whole body averaged SAR of 2-W/kg for 60 minutes of continuous RF exposure (i.e., per pulse sequence or back to back sequences/series without breaks)
<b>MR Image Artifact</b>	The presence of this implant does produce an imaging artifact. Therefore, carefully select pulse sequence parameters if the implant is located in the area of interest.

**Important Note:** The VivoKey Spark 2 is known to be magnetic (i.e., this implant has a ferrite core). However, during its intended use (i.e., implanted subcutaneously), it will not dislodge, displace, or move in association with a 3-Tesla or less, MRI environment. Importantly, the VivoKey Spark 2 becomes encapsulated by fibrous tissue within one month after implantation, thus, providing substantial stabilization of this implant.